




RIGEL

From SEAWARD





Rigel SiM Range Accessory List





rigelmedical.com
TESTED, TRUSTED... WORLDWIDE.



RIGEL SiM RANGE ACCESSORY LIST

SpO2 Interconnect Solutions (for use with SP-SiM and UNI-SiM only)				
Part No.	Description	Compatibility	Image	Specification
399A910	PULS-R - Universal SpO2 optical finger	Beijing Choice Criticare GE Tuffsat Masimo Mindray Nellcor NellcorOximax Nihon Kohden Nonin Novamatrix Phillips/HP *Custom R-curves subject to availability		Suitable with majority of SpO2 technologies
370A977	SpO2 adapter	All SpO2 adapters		Connection between SP/UNI-SiM and all SpO2 adapters








IBP Cables 1 metre length (for use with UNI-SiM only)				
Part No.	Description	Compatibility	Image	Specification
800A958	IBP Cables	Drager Siemens Axiom		10 pin male
800A959	IBP Cables	Drager Siemens Axiom		7 pin male

continues on page 3...


Please note: items not shown to scale



IBP Cables 1 metre length (for use with UNI-SiM only)				
Part No.	Description	Compatibility	Image	Specification
800A960	IBP Cables	GE Datex		10 pin male 3 key locking
800A961	IBP Cables	GE Marquette		11 pin male
800A963	IBP Cables	Mindray		5 pin male
800A965	IBP Cables	Philips Mindray		12 pin male
800A969	IBP Cables	Spacelabs, Welch Allyn		6 pin male
800A989	IBP Cables	Datascope		6 pin female
800A9101	IBP Cables	Drager		16 pin male
370A980	IBP Cables	Unterminated IBP Cables		Unterminated

Please note: items not shown to scale

Temperature Cables 1 metre length (for use with UNI-SiM only)				
Part No.	Description	Compatibility	Image	Specification
800A970	Temperature cable	Drager Siemens		7 pin male
800A971	Temperature cable	GE Datex Mindray Nihon Kohden Spacelabs Welch Allyn		Jack Ø 6.35mm
800A972	Temperature cable	GE Marquette		11 pin male
800A974	Temperature cable	Mindray		5 pin male
800A975	Temperature cable	Mindray		2 pin male
800A977	Temperature cable	Philips		2 pin male
370A981	Temperature cable	Unterminated Temperature Interconnect Cable		Unterminated

Please note: items not shown to scale





NIBP Accessories (for use with BP-SiM and UNI-SiM only)				
Part No.	Description	Compatibility	Image	Specification
370A005	Tubing sets	All manufacturers		Variety of connections and tubing to suit majority of NIBP connections

ECG Cables and Leads (for use with UNI-SiM only)				
Part No.	Description	Compatibility	Image	Specification
370A982	ECG connector box	All manufacturers		10 connections for ECG cables
298A678	Patient applied parts adapter (x11 in a pack)	All manufacturers		Connections to attached ECG connectors to ECG box

Please note: items not shown to scale

Software (for use with full SiM Range)	
Part No.	Description
383A910	Med-eBase V2.6 Asset management software
383A914	Med-eBase V2.6 additional user licence
383A916	Med-eBase V2.6 additional product licence (Rigel SiM-Range)



Printers and Scanners (for use with full SiM Range)			
Part No.	Description	Image	Specification
339A989	Bluetooth Test n Tag Elite 2 printer		Bluetooth printer for results and pass label printing (electrical safety only) (supplied with 52 x 74mm labels)
339A946	Label roll		Label roll for Test n Tag Elite 2 printer (approx. 130 off 52 x 74mm labels)
339A949	Thermal paper roll (x5 in a pack)		Thermal paper roll for Test n Tag Elite 2 printer (for results printing)
339A923	Bluetooth barcode scanner		Bluetooth scanner for asset detail transfer to Rigel products

Please note: items not shown to scale

Carry Case Solutions (for use with full SiM Range)		
Part No.	Description	Image
71G088	Grey, square carry case (suitable to hold 1 Rigel tester)	
370A983	Med-eKit Lite (Soft case suitable to hold up to 2 Rigel testers)	
Custom	Med-eKit Pro (Pelicase suitable to hold different Rigel testers, custom to requirements)	
Custom	Med-eKit Elite (Pelicase suitable to hold up to 2 Rigel testers)	
Custom	Med-eKit Plus (Durable backpack suitable to hold up to 5 different Rigel testers)	

At Rigel, we are always developing our range of accessories. If you require an accessory which is not on the list, please contact us.



Seaward, 15-18 Bracken Hill, South West Industrial Estate,
Peterlee, County Durham, SR8 2SW. United Kingdom

T: +44 (0) 191 586 3511 **E:** sales@rigelmedical.com



Rigel Medical



Rigel Medical



Rigel Medical



Rigel Medical